

1 **Title:** Evaluation of the ResistancePlus[®] MG FleXible cartridge for near point-of-care testing
2 of *Mycoplasma genitalium* and associated macrolide resistance mutations

3 **Authors:** Sweeney EL^{1#}, Mhango LP¹, Ebeyan S², Tan LY², Bletchly C³, Nimmo GR^{3,4}, and
4 Whiley DM^{1,3#}

5 ¹ University of Queensland Centre for Clinical Research (UQ-CCR), Faculty of Medicine,
6 The University of Queensland, Brisbane, Queensland, Australia

7 ² SpeedX Pty Ltd., Sydney, New South Wales, Australia

8 ³ Pathology Queensland Central Laboratory, Queensland, Australia

9 ⁴ Griffith University, School of Medicine, Queensland, Australia

10 **Running title:** Evaluation of ResistancePlus[®] MG FleXible

11 **#Address correspondence to:** Emma L. Sweeney, e.l.sweeney@uq.edu.au or David M.
12 Whiley d.whiley@uq.edu.au

13 **Keywords:** *Mycoplasma genitalium*; SpeedX; near point-of-care; macrolide; azithromycin;
14 antimicrobial resistance

15

16 **Text.**

17 *Mycoplasma genitalium* is an important sexually-transmitted infection that can cause acute
18 and/or chronic urethritis in males and females, and its prevalence is approximately 16% in
19 females and 17% in males (1). Globally, macrolide resistance is estimated to exceed 50% in
20 most urban centres (1-6), and in our local region of Queensland, Australia, exceeds 60% (7).
21 PCR is routinely used for diagnosis of *M. genitalium* and macrolide resistance mutations at
22 positions 2058 and 2059 of the 23S rRNA gene, however, the turnaround times of laboratory-

23 based methods often exceed 24hrs, and so may not be sufficiently timely to inform clinical
24 management for symptomatic patients. Here, we evaluated the performance of the
25 ResistancePlus[®] MG FleXible cartridge test on the Cepheid GeneXpert (hereafter MG-Flex),
26 which offers the potential for near point-of-care testing.

27 A bank of 181 clinical samples (*M. genitalium* positive, n = 63; and negative, n = 118) from
28 145 males and 36 females were used in this study. Original patient samples were stored at 4
29 °C for 4 weeks, then at – 20 °C for longer-term storage. Nucleic acid samples were stored at -
30 20 °C until required. Full details of the samples are provided in Supplementary Table 1. The
31 samples had all been submitted for routine *M. genitalium* testing to Pathology Queensland
32 (Brisbane, Australia) where they were tested using an in-house PCR that detects the *MgPa*
33 gene of *M. genitalium* (hereafter in-house-MgPa-PCR) (8). For the purposes of this
34 evaluation, all samples were tested using the MG-Flex near point-of-care assay (detects the
35 macrolide resistance mutations A2058T, A2058C, A2058G and A2059G) as well as the
36 Speedx ResistancePlus[®] MG test (hereafter RPMG, which detects the macrolide resistance
37 mutations A2058T, A2058C, A2058T, A2059C and A2059G). The RPMG test is a
38 Therapeutic Goods Administration (TGA, Australia) cleared and CE-IVD-marked laboratory-
39 based test.

40 The MG-Flex assay was prepared as per the kit instructions. Briefly, 44 µL of MG-Flex
41 mastermix, 1 mL of neat (unextracted) clinical sample and 10 µL of internal control were
42 added to the FleXible cartridge (Cepheid). In some instances, swab specimens (n = 17;
43 Supplementary Table 1) had less than 1 mL of sample, so the volume was made up to 1 mL
44 with sterile molecular-grade water. The RPMG assay was performed as per manufacturer's
45 instructions using stored DNA extracts from the routine testing in the in-house-MgPa-PCR.

46 The results are summarised in Table 1 and further detailed in Supplementary Table 1. The
47 MG-Flex assay detected *M. genitalium* in 61 of 63 known positive samples, providing a
48 sensitivity of 96.8%. The two samples providing negative results by the MG-Flex assay
49 (Samples 62 & 63, Supplementary Table 1) were both from males. For sample 62, a reduced
50 starting volume of sample likely contributed to discordance and in both samples, low *M.*
51 *genitalium* load was observed. 116 *M. genitalium*-negative samples were negative in the MG-
52 Flex assay, while two samples provided invalid results (Samples 180 & 181, Supplementary
53 Table 1) and there was insufficient sample to repeat testing. MG-Flex assay specificity for
54 evaluable samples was therefore 100%. The detection of macrolide resistance mutations by
55 the MG-Flex assay correlated 100% with that of the RPMG assay, with the exception of two
56 samples providing negative results for the *M. genitalium* in the MG-Flex assay (Samples 62
57 & 63, Supplementary Table 1). The overall agreement between the RPMG and MG-Flex
58 assays was 98.9%, with a kappa value of 0.98. The detection limit was also assessed by
59 testing ten-fold dilutions of a *M. genitalium* positive sample and the MG-Flex test was able to
60 reliably detect an additional dilution over the RPMG assay (supplementary Table 1).

61 In summary, we found that the MG-Flex cartridge test was highly sensitive and specific for
62 the detection of *M. genitalium* and 23S rRNA mutations.

63 **Acknowledgements:** We thank Pathology Queensland staff for their assistance with this
64 study. This study was approved by the Children's Health Queensland human research ethics
65 committee (HREC/12/QRCH/139).

66 **Conflict of Interest:** D.M.W reports research funding from Speedx Pty Ltd. The RPMG and
67 MG-Flex assays were provided by Speedx Pty Ltd. Co-authors from Speedx Pty Ltd. (SE
68 and LYT) participated in the study conception and design, data analysis and manuscript
69 writing.

70 **Funding Statement:** This work was supported by The University of Queensland strategic
71 funding.

72 **References:**

- 73 1. Getman D, Jiang A, O'Donnell M, Cohen S. 2016. Mycoplasma genitalium Prevalence,
74 Coinfection, and Macrolide Antibiotic Resistance Frequency in a Multicenter Clinical Study
75 Cohort in the United States. *J Clin Microbiol* 54:2278-83.
- 76 2. Dumke R, Thurmer A, Jacobs E. 2016. Emergence of Mycoplasma genitalium strains showing
77 mutations associated with macrolide and fluoroquinolone resistance in the region Dresden,
78 Germany. *Diagn Microbiol Infect Dis* 86:221-3.
- 79 3. Gesink D, Racey CS, Seah C, Zittermann S, Mitterni L, Juzkiw J, Jamieson H, Greer J, Singh S,
80 Jensen JS, Allen V. 2016. Mycoplasma genitalium in Toronto, Ont: Estimates of prevalence
81 and macrolide resistance. *Can Fam Physician* 62:e96-101.
- 82 4. Jensen JS, Bradshaw CS, Tabrizi SN, Fairley CK, Hamasuna R. 2008. Azithromycin treatment
83 failure in Mycoplasma genitalium-positive patients with nongonococcal urethritis is
84 associated with induced macrolide resistance. *Clin Infect Dis* 47:1546-53.
- 85 5. Gratrix J, Plitt S, Turnbull L, Smyczek P, Brandley J, Scarrott R, Naidu P, Parker P, Blore B, Bull
86 A, Shokoples S, Bertholet L, Martin I, Chernesky M, Read R, Singh A. 2017. Prevalence and
87 antibiotic resistance of Mycoplasma genitalium among STI clinic attendees in Western
88 Canada: a cross-sectional analysis. *BMJ Open* 7:e016300.
- 89 6. Deguchi T, Yasuda M, Horie K, Seike K, Kikuchi M, Mizutani K, Tsuchiya T, Yokoi S, Nakano M,
90 Hoshina S. 2015. Drug resistance-associated mutations in Mycoplasma genitalium in female
91 sex workers, Japan. *Emerg Infect Dis* 21:1062-4.
- 92 7. Sweeney EL, Trembizki E, Bletchly C, Bradshaw CS, Menon A, Francis F, Langton-Lockton J,
93 Nimmo GR, Whiley DM. 2019. Levels of Mycoplasma genitalium Antimicrobial Resistance
94 Differ by Both Region and Gender in the State of Queensland, Australia: Implications for
95 Treatment Guidelines. *J Clin Microbiol* 57.
- 96 8. Jensen JS, Bjornelius E, Dohn B, Lidbrink P. 2004. Use of TaqMan 5' nuclease real-time PCR
97 for quantitative detection of Mycoplasma genitalium DNA in males with and without
98 urethritis who were attendees at a sexually transmitted disease clinic. *J Clin Microbiol*
99 42:683-92.

100

101

102

Table 1. Summary of results across RPMG and MG-Flex assays

No. of samples	Gender (M/F)	Sample type (number)	In-house-MgPa-PCR	RPMG assay	MG-Flex assay
7	F	Cervical (3), Vaginal (1), Urine (3)	DETECTED	MG DETECTED. 23S rRNA mutations DETECTED.	MG DETECTED. 23S rRNA mutations DETECTED.
36	M	Urine (31), Rectal (3), Urethral (2)	DETECTED	MG DETECTED. 23S rRNA mutations DETECTED.	MG DETECTED. 23S rRNA mutations DETECTED.
5	F	Cervical (2), Vaginal (2), Urine (1)	DETECTED	MG DETECTED. 23S rRNA mutations not detected.	MG DETECTED. 23S rRNA mutations not detected.
13	M	Urine (11), Rectal Swab (2)	DETECTED	MG DETECTED. 23S rRNA mutations not detected.	MG DETECTED. 23S rRNA mutations not detected.
2	M	Urine (1), Rectal Swab (1)	DETECTED	MG DETECTED. 23S rRNA mutations not detected.	Mg not detected. 23S rRNA mutations not detected
92	M	Urine (85), Rectal Swab (6) Genital Swab (1)	nd	Mg not detected. 23S rRNA mutations not detected.	Mg not detected. 23S rRNA mutations not detected
24	F	Cervical (4), Vaginal (2), Urine (16), Rectal (2)	nd	Mg not detected. 23S rRNA mutations not detected.	Mg not detected. 23S rRNA mutations not detected.
2	M	Rectal Swab (1), Urine (1)	nd	Mg not detected 23S rRNA mutations not detected.	Invalid

103 nd = not detected